
**ENG / UA / PL**


**Breathalyser**

**Alco Control 8.0**

**medica+**


JAPAN TECHNOLOGY

medica+

JAPAN TECHNOLOGY

**ENG**

**medica+**  
JAPAN TECHNOLOGY

# **Instruction manual**

**Professional breathalyzer  
Alco Control 8.0**

# Contents

---

Description .....	01
Name of parts.....	02-03
Operating instructions .....	04-12
Low battery detection .....	13
Technical specifications .....	14
Packing accessories .....	15
Common faults and solution.....	15-16
Unit comparison table .....	17
Cautions.....	18
Important notes .....	19
Warranty terms.....	20-21

Dear Customer,

Thank you for choosing Medica+ products.


Please read this instruction manual carefully before using the device.

## **Humanized safe high-tech design**

This device complies with the European standard EN 16280, which sets the requirements for breath alcohol measurement devices. This standard is approved in the EU and applies to devices intended for use by the general public.

- The device features a precise and sensitive British fuel cell sensor, ensuring high accuracy, reliability, and a long service life.
- The device has a breath flow sensor that monitors the continuity of exhalation. This helps to obtain the most accurate results, close to the actual blood alcohol concentration (BAC) level.
- You can easily switch between three measurement units: mg/L, ‰BAC (per mille), %BAC (percentage in blood).

# Name of parts

## **Name of parts**

LCD – Display for information output

M Key – Unit switching and record viewing  
button

On\Off key

Mouthpiece – breath test inlet

Battery compartment

Buzzer

Vent outlet


# Operating instructions

Rinse your mouth with water before use. Avoid smoking and consuming alcohol for at least 20 minutes before using the device.

Open the battery compartment on the back of the device and insert three 1.5V (AAA) alkaline batteries.

## Power ON

After inserting the batteries, press and hold the power button for 2 seconds to turn the device on. It will automatically enter the warming-up mode.

## Power OFF


To turn off the device, press and hold the power button for 2 seconds until the "OFF", message appears. After that, the device will turn off.

In result or record viewing mode, the device will automatically turn off if not used for 30 seconds.

## Warm-up mode

Turn on the device, and it will automatically enter the warm-up state. The following image will appear to reset the sensor to ensure accurate readings.


**Note:** The warm-up time of the device varies depending on the situation:

- If the device has not been used for a long time, the warm-up will take longer.
- If the alcohol level during the last test was high ( $\geq 0.200$  mg/L), the warm-up will be longer to allow the sensor to recover and provide a more accurate result. The exact duration depends on the previous result.
- If the result was low ( $\leq 0.200$  mg/L), the device will warm up more quickly.

## Heating mode



After the warm-up, the device switches to heating mode to activate the sensor. A countdown from 19 to 0 seconds (19-00) will appear on the screen. This is necessary for the sensor to work better and become more sensitive to alcohol content.

Countdown interface showing "19" and "00"

## Breath sampling




After 19 seconds of heating, the screen will display the word "BLOW," which will flash for approximately 25 seconds along with sound signals. This indicates that the device is ready for testing. At this point, begin the measurement.

## Ready to blow interface

### How to exhale properly


Exhale steadily and continuously through the mouthpiece for about 5 seconds. A long sound signal will be heard during the exhalation. Do not stop exhaling until you see that the progress bar is filled, and the sound signal stops. After that, two short "di-di" signals will sound, indicating that the test is complete, and the result will appear on the screen.

## Dynamic Process Indicator

1. During the test, exhale slowly and evenly into the mouthpiece, without interrupting your breath, until the measurement is complete. This is crucial for obtaining an accurate result.

If you interrupt your exhale, the screen will display "Err," and you will need to repeat the test.

Err (blowing interrupted)

Test result display

## Displaying test result


If the test result exceeds the set alarm threshold, the alarm light indicator will light up red and flash 10 times, accompanied by a sound signal "di-di."

If the result is below the normal range, the screen will glow blue.

The device displays values in the range of 0.00% – 2.00% BAC:

- If the result is higher than 2.00% BAC, the screen will show "2.00% BAC."
- If the result is lower than 0.10% BAC, it will show "0.00% BAC."

With a warning level set at 0.50% BAC, the test result will be displayed as shown in the following images:


**< 0.50‰BAC**

**< 0.75‰BAC**

**= 0,50‰BAC**





**= 0,75‰BAC**

## Viewing test records

In the power-on or test result view mode, press and hold the button for 2.5 seconds to enter memory view mode. Pressing this button allows you to sequentially view up to five saved test results. The display will alternately show the record number and the corresponding test result. For example, if the second record has a value of 0.082% BAC, the screen will alternately display "-02-" and "0.082%".


The results are saved in a manner where the first record is automatically deleted when the sixth record appears, meaning the newest result overwrites the oldest.

All results are deleted from memory if the batteries are removed or the device is reset to factory settings. If no action is taken within 30 seconds in record view mode, the device will automatically turn off.

## Low battery detection

The device is powered by three 1.5-volt alkaline AAA batteries. If the battery level is low, the display will flash the message "Lob" along with the battery symbol, and a sound signal will be heard. After this, the device will automatically turn off. The batteries need to be replaced.

# Technical Specification

**Weight** – 75 g

**Dimensions** – 145 \* 58 \* 33.5 mm

**Measurement Range** – 0 ... 1 mg/L (0 ... 0.2% BAC, 0 ... 1 g/L, 0 ... 2% BAC)

**Operating Temperature Range** – 0 ... 50 °C

**Data Displayed in Three Units** – Promille, % BAC, mg/L

**Display** – Approximately 2.0-inch LCD display (blue film)

**Power Supply** – 3 AAA batteries

**Memory** – 5 test records will be stored

**Alarm Trigger Level** – 0.24 mg/L (0.05% BAC, 0.5 g/L, 0.5% BAC)

**Accuracy** – ± 5% of full scale

**Response Time** – 3-5 seconds

**Warm-up Time** – 5 seconds

**Sensitivity** – 0.001% BAC

# Packing accessories

Alcohol Tester Medica+ Alco Control 8.0 – 1 unit

Mouthpieces – 10 units

User Manual – 1 unit

**Attention!** The device does not include batteries; they must be purchased separately – 3 \* AAA alkaline batteries.

## Common faults and solutions

### Can't power on

1. Low battery – install new batteries.
2. Batteries inserted incorrectly – check the polarity and insert the batteries correctly.

### Err

1. You interrupted your exhale during the test – try again and exhale continuously for 5 seconds.

The results of two tests vary significantly in a short period.

2. You took the test too quickly after consuming alcohol – for a more accurate result, do not drink or smoke at least 20 minutes before the test.

3. Alcohol residues remain in the mouth – before the test, rinse your mouth with water or a non-alcoholic beverage to ensure the reading is closer to your actual blood alcohol level.

# Unit comparison table

No.	Blood Alcohol Content (BAC)		Breath Alcohol Content (BrAC)
	%BAC	‰BAC	mg/L
1	0.000	0.00	0.000
2	0.010	0.10	0.050
3	0.020	0.20	0.100
4	0.030	0.30	0.150
5	0.040	0.40	0.200
6	0.050	0.50	0.250
7	0.060	0.60	0.300
8	0.070	0.70	0.350
9	0.080	0.80	0.400
10	0.090	0.90	0.450
11	0.100	1.00	0.500
12	0.110	1.10	0.550
13	0.120	1.20	0.600
14	0.130	1.30	0.650
15	0.140	1.40	0.700
16	0.150	1.50	0.750
17	0.160	1.60	0.800
18	0.170	1.70	0.850
19	0.180	1.80	0.900
20	0.190	1.90	0.950
21	0.200	2.00	1.000
...	...	...	...
41	0.400	4.00	2.000

# Cautions

1. Alcohol affects everyone differently. This depends on heredity, body type, gender, and weight. Therefore, two people may drink the same amount, but their test results may differ.
2. After consumption, alcohol enters the bloodstream and gradually breaks down. As a result, the alcohol level can change within a few minutes.
3. Store the device in the box it came in. Do not leave it in direct sunlight or in damp places – this may damage it.

## Important notes

1. Do not pour alcohol or other liquids into the device – this can permanently damage it.
2. Frequent use and high alcohol levels during testing may reduce the device's sensitivity and shorten the sensor's lifespan. The more often you test high alcohol concentrations, the faster the sensor will wear out.
3. The test results are for reference only. Do not rely on them 100% to decide whether you are fit to drive. Such results may not be accepted in court.
4. For your safety and health – do not consume alcohol before driving.

The seller and manufacturer are not responsible for how you use this device.

**Remember:** do not consume alcohol before getting behind the wheel!

# Warranty Terms

Dear customers,

Please keep this instruction manual properly to use the warranty card. The warranty is provided for one year from the date of purchase (or for 2 or 3 years if the extended warranty is purchased) in case of malfunction under normal operating conditions.


**The warranty does not cover the following cases:**

1. If the product is damaged due to improper use or lack of proper maintenance.
2. If the product is damaged as a result of disassembly or modifications.
3. If parts not supplied by our company were applied to and used with the device.
4. Damage caused by force majeure circumstances.

According to the terms, technical maintenance provided outside the warranty is subject to additional charges.

Parts of the product replaced during the warranty period remain the property of MEDICA+.

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