ENG / UA

Ultrasonic Toothbrush

PROBRUSH 7.0

medica+
JAPAN TECHNOLOGY

medica+

JAPAN TECHNOLOGY

ENG

Instruction Ultrasonic Toothbrush

Probrush 7.0

Content

What's Included	01
Features	02
Safety first.....	03-04
Contents	05
Before Use	06
Brush Head	07
Travel Cover/Mount	08-09
Brushing Instructions	10
Quadrant Method	11
Care Information	12
Replacement Brush Heads.....	13
Troubleshooting.....	14

Dear customer,
Thank you for choosing MEDICA+
products.

Please read this instruction before use
and use it in a proper way so as to
achieve ideal effects

What's Included

- Tooth brush handle ×1
- Brush head ×2
- Brush head cap ×1

Features


1. Cleans effectively with 24,000 vpm.
2. Ergonomic design, non-slip, and comfortable to hold.
3. Waterproof Performance: IPX7
4. Battery Life: Up to 120 days (when used twice per day at 2 minutes each time).
5. Notification every 30 seconds to switch to next quadrant of mouth.
6. Auto-Shutoff after 2 minutes
7. 2 modes: Clean / Massage
8. Travel Cover/Mount: Attach to mirror or a wall or use it to protect the brush head when traveling

Safety first

1. Do not open, repair, or alter the product without authorization.
2. Do not knock or drop the product.
3. If product is defective, please stop using and see the Troubleshooting section in this manual.
4. Only use with MEDICA+ brush heads.
5. If you've had oral surgery during the past two months, please consult with your dentist before using this toothbrush.
6. If your gums bleed after using this toothbrush and bleeding continues for more than one week, please stop use and consult your dentist before beginning use again.
7. This product is designed to only clean your teeth and tongue. Do not use for any other purpose.



8. Children should be supervised to ensure proper usage.
9. Keep the toothbrush away from high temperatures and direct sunlight.

Contents


Before Use

1. Hold the brush handle with the battery cover facing toward you.
2. Loosen the battery cover by turning it counterclockwise.
3. Pull the battery cover away from the brush handle.
4. Remove the clear plastic tab and close the battery cover tightly.


Brush Head

Pull the brush head away from the handle to remove it or push it down over the metal stem to attach it to the brush handle.


Travel Cover/Mount

1. Peel the backing paper from the adhesive strip.

2. Mount to a mirror or wall by pressing the sticker side of the cover firmly against the surface of a smooth clean wall or mirror.*


* Sticker may not adhere to all surfaces. Make sure surface is clean and fully dry beforehand.

Brushing Instructions

1. Wet the bristles and apply an appropriate amount of toothpaste or whitening foam.
2. Place the toothbrush bristles against the teeth.
3. Press the button and choose your mode:
 - 1st press: Clean Mode
 - 2nd press: Massage
 - 3rd press: Off
4. Gently move the brush head slowly across the teeth in a small up and down motion. For best results, use the Quadrant Method, brushing each quadrant for 30 seconds (there will be a pause in vibration when it's time to switch).

Quadrant Method

Auto shut-off occurs after 2 minutes.

When finished, rinse the toothbrush and put it in a ventilated and dry place.

Care Information

1. Rinse brush head off thoroughly with water after every use. Use a cloth to wipe of the handle and protective cover. A build up of debris could cause your toothbrush to not function properly and this is not covered by the warranty.
2. Replace the battery when low power is observed. Do not leave a battery in the handle past its life. Battery corrosion is not covered by the warranty.

Replacement Brush Heads

To be used only with MEDICA+ replacement brush heads.

We recommend replacing the brush heads every 2-3 months.

Troubleshooting

Case	Possible Causes	Solution
Toothbrush has a weak vibration.	• Too much force on the bristles against teeth. • Low battery.	• Use less pressure when brushing. • Replace battery.
Toothbrush loses power before 2 minutes.	• Low battery.	• Replace battery.
Toothbrush does not power on.	• Low battery. • Other cause.	• Replace battery. • Contact us.

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