ENG / UA/ PL

ICE PRO 10X

medica+ 
JAPAN TECHNOLOGY

ENG

Instruction

MEDICA+ ICE PRO 10X
cooling photoepilator

Content

Safety Instructions	01-02
Device Description	03-07
Contraindications	08-09
Expected Results	10
Preparing the Device for Use.....	11
Instructions for Use.....	11-17
Replacing the Attachment.....	18
Error Notification Function.....	19
Use of Attachments	19-20
Cleaning the Device	21
Warranty Terms.....	22
Technical Specifications	23

Dear Customer,
Thank you for choosing MEDICA+ products.
Please read this instruction manual carefully
before using the device.


Safety Instructions

1. The device lamp does not support connection or disconnection while energized. Before replacing the attachment, make sure to switch off the device and disconnect it from the power supply. If the attachment is not installed, the device will not turn on. During use, ensure that the ventilation openings remain unobstructed.
2. Do not use the device on the face above the cheek line, particularly in the area of the eyes, eyebrows, and eyelashes, as this may cause serious eye injury.
3. Since the device operates from an external power source, it is prohibited to use it at a distance far from the electrical outlet or with a stretched power cable.
4. The active treatment area is the light emission window. During operation, the temperature in this area may reach up to 45°C. To control temperature exposure, it is recommended to limit contact with one area of skin to no more than 1 minute per session.
5. Protective goggles must be worn when operating the device.
6. Before the procedure, cleanse the area of skin to be treated. Do not apply cosmetic or skincare products during use of the device in order to prevent damage to the device and reduced effectiveness of the procedure.


7. During first use, start at the lowest intensity level and gradually select the optimal level without causing harm to the skin.
8. If redness or swelling appears after using a high intensity level, it is recommended to use the cold compress function or additionally cool the treated area.
9. After completing the procedure, it is recommended to moisturize the skin and ensure reliable sun protection, avoiding exposure to intense light and ultraviolet radiation.
10. After switching off the device and disconnecting it from the power supply, carefully wipe the lamp connector with a damp cloth or soft tissue. Avoid allowing water or other liquids to enter the device, as this may cause damage.
11. Do not allow the device or its accessories to come into contact with liquids and do not use them in a humid environment.
12. Store the device in a dry, well-ventilated place out of the reach of children.
13. Do not attempt to disassemble or repair the device yourself.
14. If the device does not function properly, deliver it together with its accessories to a MEDICA+ service center for maintenance by a qualified specialist.

Device Description


The MEDICA+ ICE PRO 10X is a next-generation high-tech IPL photoepilator with an active ice-cooling system, designed for effective, fast, and максимально comfortable removal of unwanted hair at home.


Flash button

Power Adapter

1. Treatment Window
2. Skin Tone Sensor
3. Electronic Display
4. Indicator Light
5. Power / Intensity Adjustment Button

1. Flash Attachment Type Icon
2. Automatic Mode Icon
3. Cooling Mode Icon
4. Flash Count Indicator
5. Intensity Level Indicator

1

2

3

4

1. Manual Mode
2. Automatic Mode
3. Automatic Skin Tone Detection Mode
4. Skin Tone Detection Mode

Power Button (ON/OFF) and Intensity Level Selection

The button is used to turn the device on and off, as well as to select the treatment intensity level.

1. Press and hold the button for more than 2 seconds (after connecting the device to the power source) to turn it on. The device starts at Level 1.
2. Each short press of the button after powering on changes the intensity level — from lower to higher. A total of 5 treatment levels are available.
3. To turn off the device, press and hold the power button again for more than 2 seconds.

Ice Cooling Function

1. The device is equipped with an active ice cooling function that automatically activates when the device is powered on.
2. Press and hold the Flash Button to activate or deactivate the ice cooling function.
3. After deactivation, the corresponding cooling icon on the display will turn off.

Note:

The ice cooling function is intended to reduce skin temperature during or after the hair removal procedure, minimize the sensation of heat, and enhance overall comfort. It does not have any medical or therapeutic effect.

Flash Button

The flash button activates a light pulse only when the treatment window is in full contact with the skin.

1. A short press triggers a single light pulse.
2. A long press is used to activate or deactivate the ice cooling function.

Button Combination

1. Flash Button + Power Button (short press):
Switches between Normal Mode, Automatic Flash Mode, and Automatic Skin Tone Detection Mode.

Indicator Status	Mode	Meaning
Off	Off	The device is not operating (not connected to a power source)
Slow blue flashing	Activation / Standby	The device is powered on and in standby mode
Fast blue flashing	Level Preparation	The device is preparing to operate at the selected intensity level
Steady blue light	Operating Mode	The device is ready to deliver a pulse when full skin contact is achieved

Contraindications

1. The device must not be used by individuals under the age of 16.
2. Use is contraindicated for pregnant women, women planning pregnancy, or during breastfeeding due to the increased risk of skin irritation or damage.
3. Do not use the device if you have cardiovascular diseases, high blood pressure, malignant tumors, infectious skin diseases, burns, or blood clotting disorders.
4. Do not use the device on areas with sunburn or other skin injuries until fully healed. Use on damaged skin increases the risk of burns.
5. The device is contraindicated for individuals with epilepsy, as bright light flashes may worsen the condition.
6. Do not use the device on the head, ears, nipples, or genital area, as this significantly increases the risk of burns and skin injury.
7. Do not use the device on areas with tattoos, dark brown or black spots, including dark freckles, moles, birthmarks, warts, etc. Such formations may absorb excessive light energy, increasing the risk of skin burns.
8. Use is prohibited in the presence of pain, psoriasis, eczema, burns, herpes, open wounds, or rashes in the treatment area. In such cases, application may lead to increased skin sensitivity, swelling, blistering, or infection. Consult a physician before use and ensure that the skin has completely healed.

9. The device must not be used by individuals prone to scarring, with photosensitivity, or by those taking medications that increase skin sensitivity to light, including but not limited to: nonsteroidal anti-inflammatory drugs (aspirin, ibuprofen, paracetamol), tetracyclines, phenothiazines, diuretics, sulfonamides, triazines, imidazolamide, fluorouracil, vinblastine, griseofulvin, retinoids (including both topical and systemic forms), carotenoids, etc.
10. Do not use the device if you have photosensitive conditions such as porphyria, polymorphic light eruption, or solar urticaria, as these conditions increase skin sensitivity to light.
11. Use of the device is not recommended for women during menstruation.
12. Do not use the device if you have implanted electronic devices such as a pacemaker, insulin pump, urinary incontinence control device, etc.
13. Use is prohibited in cases of frequent herpes recurrences in the treatment area.

Expected Results

Immediately After the Procedure

After treatment, significant side effects are generally not observed.

During the First Few Weeks After Starting Treatment

In the first weeks following the initial sessions, you may notice continued growth of some hairs. These are hairs that were in the resting phase (not in the active growth phase) at the time of treatment.

After 6 Weeks Within the 12-Week Course

A noticeable reduction in hair growth can be observed. However, some hairs may not yet have been treated during their active growth phase, so it is important to continue regular weekly treatments.

After Completing the 12-Week Course

After completing the full treatment course, a significant reduction in hair in the treated area is typically observed. The remaining hair becomes finer and lighter. To maintain the results, it is recommended to perform maintenance treatments once a month or as needed.

Preparing the Device for Use

1. Remove all components of the device from the packaging.
2. Before each session, inspect the attachment window and the power cable for any damage.
3. Before turning on the device, ensure that the light attachment is properly and securely installed.
4. Connect the power cable to an electrical outlet and the other end to the corresponding port on the device. After connection, the indicator light at the top of the device will slowly flash blue, indicating standby mode.

Instructions for Use

Determine your skin type according to the phototype scale to ensure the device can be used safely. Before the procedure, make sure the skin is healthy, without damage, irritation, or open wounds. The treatment area must be properly prepared, clean, and dry, free from cosmetic products, creams, makeup, deodorants, or antiperspirants.

Skin color Hair color	White	Ivory	Beige	Light Brown	Brown	Dark Brown
 Black						
 Dark brown						
 Brown						
 Golden						
 Red						
 White						
 Gray						

 : Indicates applicable :  Indicates not applicable

If you are using the device for the first time or plan to treat a new area of the body, be sure to perform a skin sensitivity test.

1. Select a small area of skin close to the intended treatment area.
2. Set the intensity to the lowest level.
3. Following the instructions for use, deliver only 2-3 light pulses.
4. After the test, wait 24 hours and assess the condition of the skin.

If after 24 hours the skin appears normal, without redness, irritation, or other adverse reactions, you may proceed with the full treatment.

Recommended Treatment Protocol

The hair growth cycle consists of three phases: anagen (active growth phase), catagen (transitional phase), and telogen (resting phase). At any given time, not all hairs are in the same phase. IPL technology is most effective during the anagen phase — when the hair is actively growing and connected to the hair follicle. Therefore, it is recommended to perform treatments regularly over a 12-week period to ensure that each hair is treated at least once during its growth phase.

The first 4 treatments are performed once a week. This allows for maximum effectiveness on hairs in the active growth (anagen) phase.


From the 5th to the 12th treatment, sessions are performed once every 2 weeks. This interval allows treatment of hairs that enter the growth phase later.

After completing the 12-treatment course, to maintain results it is recommended to perform 1 session every 1-2 months or as needed if occasional hairs appear.


Before the Procedure

Before the procedure, completely shave the hair in the treatment area using a razor, preferably one day prior to the session. The skin should be clean and smooth. Do not use wax, an epilator, or depilatory creams, as they remove hair from the root and reduce the effectiveness of IPL. Before first use or when treating a new area, perform a skin sensitivity test.


1. Press and hold the ON/OFF button for 2 seconds to turn on the device. The fan will start operating, and the electronic display will light up. To activate the ice cooling function, press and hold the Flash button for 3 seconds. The cooling icon will appear on the screen, and the area around the light emission window will become cold.

2. The device has three operating modes: Manual Mode, Automatic Mode, and Automatic Skin Tone Detection Mode. In Manual Mode, the user presses the Flash button for each pulse. In Automatic Mode and Automatic Skin Tone Detection Mode, the device automatically emits light pulses after contact with the skin. To switch between modes, briefly press the On/Off + Flash buttons simultaneously.


3. Hold the photoepilator perpendicular to the skin at a 90° angle. When the attachment is in full contact with the skin, the indicator will turn blue and remain steadily lit. This indicates that the device has correctly detected the skin and is ready for operation.


4. After switching on, the device automatically sets the intensity level to Level 1. By briefly pressing the On/Off button, you can adjust the level from 1 to 5 according to your skin tolerance.

Note: Each time a flash is emitted, you may feel slight warmth or a tingling sensation — this is a normal skin reaction during the procedure.

5. To remove hair, press the Flash button. After each pulse, move the attachment to an adjacent area of skin. Allow pauses for the skin to cool, and if necessary, repeat the treatment 1-2 times on the same area. Do not treat the same area repeatedly without allowing it to rest in order to avoid burns. Do not perform more than three consecutive passes on the same area during a single session.

6. In Automatic Skin Tone Detection Mode, the sensor determines the skin tone and automatically selects the appropriate intensity level. Ensure that the sensor is in full contact with the skin. If contact is incomplete, the energy level will be 0 and no pulses will be emitted.
7. After hair removal, press and hold the ON/OFF button for more than 2 seconds to turn off the device. The fan will stop operating after 10 seconds, after which the device can be disconnected from the power supply.

Note: If the device is not used for 5 minutes, it will automatically enter standby mode: the indicator will slowly flash blue, and the fan will stop.

8. After the device has cooled down, clean the treatment area. Then place the device and attachments back into the box and store them in a dry, cool place.
9. Skin Soothing: If necessary, apply a soothing product, such as MEDICA+ Fluid After Epilation, to reduce possible temporary redness or warmth. Regularly moisturize the skin for several days after the procedure and avoid direct sun exposure. Use sunscreen if necessary.


Replacing the Attachment

The standard lamp attachment is designed for 500,000 light pulses. When the remaining pulse counter displays "000000," it means that the lamp attachment has reached the end of its service life and must be replaced.

To purchase or replace the lamp attachment, please contact the official seller or the manufacturer.


Removing the Lamp Attachment:

Hold the lamp attachment and carefully pull it in the direction indicated in the illustration to detach it from the device housing.

Installing the Lamp Attachment:

Align the contact element of the lamp attachment with the corresponding connector on the device, insert the attachment, and press firmly until you hear a click confirming proper installation.

Error Notification Function

Er 01 — Fan Error

Er 03 — Temperature Error

Solution:

If any of the above error codes appear, turn off the device, disconnect it from the power supply, wait a few minutes, and then turn it on again. If the error persists, discontinue use of the device and contact the service center.

Use of Attachments

HR Attachment (Hair Removal) – Yellow-Tinted Lamp (Included)

The HR attachment features a wavelength range of 510–1200 nm and is designed for effective removal of unwanted hair on various areas of the body. It is suitable for use on the legs, arms, underarms, back, chest, bikini area (external area only), and on the face below the cheekbone line.

The light pulse targets melanin in the hair follicle, contributing to a gradual reduction in hair growth intensity with regular use.

Do not use the device around the eyes or directly on the genital area.

SR Attachment (Skin Rejuvenation) – Orange-Tinted Lamp (Sold Separately)

The SR attachment features a wavelength range of 560–1200 nm and is used to improve overall skin condition, enhance skin firmness, and contribute to a more even skin tone while reducing the appearance of pigmentation.

It is suitable for treatment of the face (below the cheekbone line), neck, décolleté area, and other areas requiring skin tone correction.

Avoid use around the eyes and on mucous membranes.

AC Attachment (Acne Treatment) – Clear Lamp (Sold Separately)

The AC attachment features a wavelength range of 400–700 nm and is designed to reduce inflammatory skin conditions and the appearance of acne.

It may be used on areas with breakouts, including the face (below the cheekbone line), back, and shoulders. Do not use the attachment around the eyes, on damaged skin, or on mucous membranes.

To prevent discomfort in sensitive areas such as the bikini area, inner thighs, and buttocks, it is recommended to begin the procedure at the lowest intensity level.

Cleaning the Device

1. Before cleaning, turn off the device and disconnect it from the power supply.
2. If stains or dark marks are present, wipe the treatment window with a dry, lint-free cloth.
3. Do not use acetone, alcohol-based, or abrasive products, as they may damage the device.
4. Keep the device dry. Do not immerse the device or any of its parts in water, as this may result in electric shock.
5. After cleaning, perform a visual inspection of the device. If dark marks remain, repeat the cleaning process until they are completely removed.

Maintenance and Storage

To prevent damage or loss of components, it is recommended to store the device and its accessories in the original packaging after each use.

Before storage, clean the device and place it in a safe, dry, and protected location.

Warranty Terms

Dear Customers,

Please keep this instruction manual properly in order to use the warranty certificate. The warranty is valid for one year from the date of purchase (or for 2 or 3 years if an extended warranty has been purchased) in the event of a malfunction under normal operating conditions.

The warranty does not cover the following cases:

1. If the product is damaged due to improper use or lack of proper maintenance.
2. If the product is damaged as a result of disassembly or modification.
3. If parts not supplied or approved by our company have been used with the device.
4. Damage caused by force majeure circumstances.

According to these terms, maintenance services provided outside the warranty period are subject to additional charges.

Parts replaced during the warranty period remain the property of MEDICA+.

Technical Specifications

Device Lifespan: 500,000 pulses

Light Type: Intense Pulsed Light (IPL)

Light Source: Xenon flash lamp

Wavelength Range: 510–1200 nm

Pulse Duration: 0.4–5.5 ms

Treatment Window Area: 3.9 cm²

Intensity Levels: Levels 1–5

Maximum Pulse Energy: 23 J $\pm$ 20%

Energy Density: 1.02–5.90 J/cm²

Pulse Interval: 0.5–2 s

Input Voltage: AC 100–240 V

Rated Current: 0.1–2.5 A

Rated Frequency: 50/60 Hz

Service Life: 5 years

Operating Temperature: 0°C–30°C

Operating/Storage Humidity: 25%–75%

Operating/Storage Atmospheric Pressure: 80–106 kPa

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