ENG / UA / PL / DE

ProBrush 9.0

medica+
JAPAN TECHNOLOGY

ENG

medica+
JAPAN TECHNOLOGY

Instruction

MEDICA+ PROBRUSH 9.0 Toothbrush

Contents

Safety Instructions 01-02

Device Description 03-04

Control Panel 05-06

Operating Modes 06-07

Instructions for Use 08-09

Installing the Brush Head and Operation.... 10-11

Cleaning and Storage 12

Device Charging 13

Warranty Terms 14

Technical Specifications 15

Dear Customer!

Thank you for choosing MEDICA+ products. Please read this instruction manual carefully before using the device.

Safety Instructions


1. This device is intended for individual use. The toothbrush should not be used by persons with reduced physical, sensory, or mental capabilities unless they are supervised or have received proper instructions regarding safe use.
2. Persons who have undergone oral or gum surgery within the last two months are advised to consult a dentist before starting to use the product.
3. The charger must not be used outdoors, in humid environments, or near heat sources and heated surfaces. Do not allow the charger to come into contact with water.
4. The product is intended exclusively for hygienic cleaning of teeth, gums, and tongue. Use of the product for purposes other than its intended use is prohibited. In case of pain, irritation, or discomfort, stop using the product immediately and consult a doctor or dentist.

5. A slight sensation of irritation during the first use of the device is normal. If irritation persists for more than two weeks, discontinue use and consult a dentist.
6. Do not use brush heads with damaged, deformed, or worn bristles. The recommended replacement interval for the brush head is at least once every 3 months or earlier if signs of wear appear.
7. When using toothpastes containing sodium bicarbonate, hydrogen peroxide, or other abrasive components (for example, whitening toothpastes), the brush head must be thoroughly rinsed with running water after each use.
8. Slight heating of the device housing during charging is acceptable. If the device overheats, disconnect it from the power supply and contact the MEDICA+ service center for inspection.
9. The manufacturer is not responsible for damage to the product or injury caused to the user as a result of improper use or failure to comply with the requirements of this instruction manual.


Device Description

The MEDICA+ PROBRUSH 9.0 toothbrush with high-frequency sonic vibration combines high cleaning efficiency with modern oral care technologies. With a frequency of up to 38,000 vibrations per minute, it gently yet thoroughly removes dental plaque, food debris, and bacteria even in hard-to-reach areas. The brush head features a soft rubber coating that ensures gentle and safe cleaning of sensitive teeth and gums. The special rubber coating protects enamel from possible damage during vibration, reducing stress on the teeth and increasing comfort during daily use.


The device is equipped with 5 operating modes, allowing the brushing process to be individually adjusted according to different needs — from daily care to a gentle mode for sensitive teeth and gums. Three intensity levels make it possible to select the optimal vibration strength for maximum comfort.


Main unit (toothbrush) — 1 pc.


Travel case — 1 pc.

Charging base — 1 pc.


Replacement brush head — 2 pcs.

Type-C charging cable — 1 pc.

Control Panel

Power Button

- A short press turns the toothbrush on and allows you to select the desired mode.
- To switch off the toothbrush, press and hold the button for 4 seconds.

Intensity Adjustment Button

By pressing the intensity adjustment button, you can select the vibration intensity according to your personal preferences to ensure the most comfortable and effective cleaning.

Intensity Indicator

The device is equipped with three intensity levels: low, medium, and high.

Mode Indicator

Displays the selected mode: sonic vibration, cleaning, whitening, sensitive, and massage. During brushing, the user can change the mode at any time by briefly pressing the power button to switch to the next mode.

Charging Indicator

Displays the charging status of the device: charging, fully charged, low battery level.

Operating Modes

After switching on the toothbrush, press the adjustment button to increase or decrease the vibration intensity. Select a comfortable level according to your personal preference.

Sonic Vibration: 36,000 times/min

A high-frequency sonic vibration mode designed for comprehensive tooth cleaning. Provides an optimal balance between cleaning efficiency and gentle impact on enamel and gums. Suitable for daily use and maintaining proper oral hygiene. Recommended for users who prefer gentle yet thorough cleaning.

Cleaning: 38,000 times/min

An intensive mode for deeper cleaning. Suitable for users who are already accustomed to sonic toothbrushes or require more powerful cleaning.

Whitening: 28,000–38,000 times/min

A mode with rhythmic alternation of high and low vibration frequencies. Promotes effective tooth cleaning and helps remove plaque and food debris.

Sensitive: 28,000–38,000 times/min

A mode designed for sensitive teeth and gums. The vibration intensity gradually changes from low to higher levels, ensuring smooth adaptation and gentle cleaning without irritation.


Massage: 28,000–38,000 times/min

A gum massage mode with cyclic changes in vibration intensity (from strong to weak and back). Helps improve blood circulation in the gums and provides additional oral care.

The toothbrush automatically memorizes the last selected mode and intensity level. After switching on, the device immediately resumes operation using the previously selected settings, allowing you to continue brushing with your preferred configuration.

Instructions for Use


To maintain proper oral hygiene, it is recommended to brush your teeth twice a day for 2 minutes. The device is equipped with a built-in timer that helps ensure the correct brushing time: during operation, a short pause occurs every 30 seconds, signaling that it is time to change the brushing area. This allows all areas of the teeth to be cleaned evenly and helps achieve maximum cleaning effectiveness.


Recommended Brushing Sequence

- Outer surfaces of the teeth: upper and lower jaws — left side;
- Inner surfaces of the teeth: upper and lower jaws — left side;
- Inner surfaces of the teeth: upper and lower jaws — right side;
- Outer surfaces of the teeth: upper and lower jaws — right side.


Installing the Brush Head and Operation


1. Carefully install the brush head onto the toothbrush body. Leave a small gap between the brush head and the body — this is necessary to ensure proper vibration. After installation, switch on the device and make sure it operates correctly.


2. Apply a small amount of toothpaste to the bristles. Do not switch on the device before this step to avoid splashing the toothpaste.


3. Lightly moisten the brush head and toothpaste before brushing. This helps the toothpaste spread more evenly and clean the teeth more effectively.

4. Brush your teeth using gentle movements without applying excessive pressure to avoid injuring the gums and enamel. Move the brush head slowly back and forth across the teeth.

5. Do not press the brush against the teeth during brushing; allow the toothbrush to do the work for you.

6. After completing the procedure, rinse your mouth with water to remove any remaining toothpaste and debris.

Cleaning and Storage

Regular cleaning of the toothbrush ensures safe operation and extends the service life of the device.

Cleaning the Brush Head

After each use, clean the brush head and allow it to dry naturally. Avoid applying excessive pressure or pulling the bristles to prevent deformation and shedding.


Cleaning the Main Unit

Wipe the toothbrush body with a soft cloth or a damp tissue and allow it to dry completely. Do not immerse the body in water for extended periods and do not use hot water or aggressive cleaning agents.

Cleaning the Charger

Wipe the charger with a soft cloth and allow it to air dry. Do not immerse the USB port in liquid and do not clean it with water or acidic solutions.


Replacing the Brush Head

Wipe the charger with a soft cloth and allow it to air dry. Do not immerse the USB port in liquid and do not clean it with water or acidic solutions.

Device Charging

1. Before first use, fully charge the device.
2. Connect the charging base to a power source.
3. Place the toothbrush on the charging base. Charging occurs automatically via the inductive contact.
4. During charging, the indicator lights up red. After full charging, the indicator changes to green.
5. If the device has not been used for a long time, the indicator may light up with a delay — this is normal device operation.
6. Following the charging instructions ensures the maximum service life of the battery.

Warranty Terms

Dear Customers!

Please keep this instruction manual properly in order to use the warranty card. The warranty is provided for one year from the date of purchase (or for 2 or 3 years in the case of purchasing an extended warranty) in the event of a malfunction under normal operating conditions.

The warranty does not cover the following cases:

1. If the product is damaged due to improper use or lack of proper maintenance.
2. If the product is damaged as a result of disassembly or modifications.
3. If parts not belonging to our company were applied to or used with the device.
4. Damage caused by force majeure circumstances.

According to the terms, technical service provided outside the warranty period is subject to additional payment.

Parts of the product replaced during the warranty period remain the property of MEDICA+.

Technical Specifications

Device weight: 0.135 kg
Dimensions (H × W × D): 18 × 2.7 × 2.7 cm
Housing material: ABS + PP
Water resistance rating: IPX7 (protection against moisture and short-term immersion in water)
Battery type: lithium
Battery capacity: 2000 mAh
Charging time: approximately 20 hours
Operating time: up to 30 days
Charging type: inductive (via charging base)
Input voltage: AC 100–240 V
Charger output: 5 V / 50 mA
Power: 2.5 W
Number of operating modes: 5
Brushing timer: 2 minutes
Vibration frequency: up to 38,000 vibrations/min
Number of movements: up to 80,000 movements/min
Noise level: ≤ 75 dB (A)
Bristle type: DuPont

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