ENG / UA/ PL

Vibroskin 9.0

medica+
JAPAN TECHNOLOGY

ENG

medica+
JAPAN TECHNOLOGY

Instruction manual

MEDICA+ Vibroskin 9.0
Ultrasonic Facial Skin Scrubber

Content

Safety Precautions..... 01-02

Suitable for skin with..... 02

Description..... 03-04

Contraindications..... 05-06

Recommendations for Use..... 06

Operating Modes..... 07-13

Switching On and Off..... 13

Charging the Device..... 14

Cleaning the Device..... 15

Warranty Terms 16

Technical Specifications..... 17

Dear Customer!

Thank you for choosing Medica+ products. Please read these instructions carefully before using the scrubber.

Safety Precautions

1. Fully charge the device before using it for the first time.
2. Before charging the device, make sure that the charging cable is not damaged or deformed.
3. Before storing the device, make sure that it is switched off. Store the device in a dry place out of the reach of children.
4. Use antiseptic wipes or a cotton pad moistened with chlorhexidine solution to clean the housing and the ultrasonic scrubber spatula. Do not use abrasive, aggressive or corrosive cleaning agents. Clean the device after each use.
5. Do not throw the device into a fire or store it in high-temperature or humid environments.
6. Do not immerse the device in water or rinse it directly under running water.

7. Do not use the device if it has been damaged as a result of being dropped.
8. Do not disassemble or repair this product yourself. For repairs, please contact a Medica+ service centre.

Suitable for skin with

- dull complexion;
- enlarged pores;
- excess sebum;
- uneven texture;
- tendency to develop comedones;
- reduced firmness;
- signs of photoaging.

Description

MEDICA+ Vibroskin 9.0 is a multifunctional cosmetic device for comprehensive skin care. It combines advanced cleansing and rejuvenation technologies to help maintain healthy-looking, well-cared-for skin.


In Sonic mode, the device uses ultrasonic vibrations at a frequency of 26 kHz, providing ultrasonic cleansing of the skin surface, gentle exfoliation of dead skin cells, and removal of surface impurities and excess sebum.

In addition to ultrasonic technology, the scrubber features the following functions:

- Ultrasonic cleansing;
- Lifting to tone the skin and help maintain its firmness;
- Iontophoresis to improve the penetration of active ingredients in cosmetic products;
- Desincrustation for more effective pore cleansing;

1. Scrubber spatula
2. LED indicator light
3. Ion– mode (moisturizing / nourishing)

4. Ion+ mode (deep cleansing)
5. Lifting mode (lifting)
6. Sonic mode (ultrasonic cleansing / peeling)
7. On/Off button
8. Battery charge indicator

The combination of these technologies promotes deep skin cleansing, improves microcirculation, enhances skin firmness and elasticity, stimulates natural regeneration processes, reduces imperfections, and improves the overall condition of the skin.

Contraindications

Although the ultrasonic skin scrubber is a safe cosmetic device for home use, its use is not recommended in the following cases:

- pregnancy and breastfeeding;
- presence of a pacemaker, implanted defibrillator, or other electronic medical implants;
- cardiovascular diseases or heart rhythm disorders;
- epilepsy and other neurological conditions for which the use of electrical cosmetic devices is contraindicated;
- cancer;
- acute infectious diseases or elevated body temperature;
- damaged skin, open wounds, burns, or severe inflammation in the treatment area;
- active dermatological conditions in the treatment area (dermatitis, eczema, psoriasis, etc.);
- severe allergic reactions or exacerbation of skin conditions;
- recent chemical peeling, intensive cosmetic procedures, or surgery — until the skin has fully recovered.

Caution! Do not use the device around the eyes, on the eyelids, mucous membranes, or damaged skin.

The device is not intended for use by children or by persons with physical, sensory, or cognitive impairments without supervision by a responsible person. Before first use, if you have any chronic medical conditions or concerns, it is recommended that you consult a physician.

Recommendations for Use

1. Before use, remove the protective cap from the ultrasonic spatula.
2. Press and hold the power button to switch the device on or off.
3. After switching the device on, briefly press the button (approximately 1 second) to cycle through the operating modes (Ion– / Ion+ / Lifting / Sonic).


Operating Modes

Mode 1 — Ion– (Moisturizing Mode) + Blue Light

This mode is used for gentle massage and helps improve the penetration of cosmetic products during the procedure.
High-frequency vibrations stimulate blood microcirculation and activate cellular metabolism. Regular use helps make the skin firmer, smoother, and well cared for.

Purpose: to promote the penetration of active ingredients in cosmetic products during the procedure, stimulate microcirculation, tone the skin, and improve its overall condition.
In this mode, hold the scrubber with the flat side of the spatula against the skin. Perform the procedure for 5 minutes, slowly moving the scrubber from the center of the face outward.


The procedure is performed using MEDICA+ LEDerm Blue Booster, serums, creams, essences, and face masks.
In this mode, the scrubber can be used 3–4 times a week.


Mode 2 — Ion+ (Deep Cleansing) + Red Light

This is an electrical mode used for deep cleansing / desincrustation. It helps cleanse the pores of excess sebum and impurities that accumulate within them. Cleansing is achieved through ionic/electrical action.

Purpose: deep cleansing of the pores and removal of excess sebum and impurities that accumulate within them.
In this mode, hold the scrubber with the edge of the spatula against the skin at a 45-degree angle. Move the spatula from the outer areas of the face toward the center.

Before the cleansing procedure, it is recommended to apply MEDICA+ Hydration Gel to the skin and leave it under plastic film for 15 minutes. After this time, remove the film from the area to be treated and begin the cleansing procedure. If necessary, the skin can be additionally moistened using a cotton pad soaked in toner or warm water.

After completing the procedure, apply MEDICA+ Soothing Mask or MEDICA+ Pore Control Mask to the skin.

- For oily and combination skin, the procedure can be performed 1–2 times a week.
- For dry facial skin — once a week.


The duration of the procedure is individual and depends on the degree of skin congestion. On average, it should not exceed 10–15 minutes.

Mode 3 — Lifting + Blue Light

This mode is used to complete the skincare procedure, tone the skin, and help maintain a feeling of firmness and smoothness.

Purpose: to complete a comprehensive skincare routine, tone the skin, help maintain its firmness and smoothness, and soothe the skin.
In this mode, hold the scrubber with the flat side of the spatula against the skin. Slowly move the scrubber from the center of the face outward.


The procedure is performed using MEDICA+ LEDerm Blue Booster, serums, creams, essences, and face masks.

- For an intensive course: 3–4 times a week.
- To maintain the results: 1–2 times a week.

Mode 4 — Sonic (Ultrasonic Cleansing) + Red Light

Provides ultrasonic cleansing of the skin surface. Ultrasonic vibrations at a frequency of 26 kHz gently remove dead skin cells. The procedure refreshes the skin, evens out its texture, improves the complexion, and enhances the effectiveness of subsequent skincare. Ultrasonic vibrations gently remove dead skin cells, surface impurities, makeup residue, and excess sebum. Helps make the skin smoother, fresher, and more even in texture.

Purpose: surface cleansing, gentle exfoliation of dead skin cells, and improvement of skin texture. Dead skin cells → surface impurities → smoothness → more even texture.
In this mode, hold the scrubber with the edge of the spatula against the skin at a 45-degree angle. Move the spatula from the outer areas of the face toward the center.

Before the cleansing procedure, it is recommended to apply MEDICA+ Hydration Gel to the skin and leave it under plastic film for 15 minutes. After this time, remove the film from the area to be treated and begin the cleansing procedure. After completing the procedure, apply MEDICA+ Soothing Mask or MEDICA+ Pore Control Mask to the skin.

- For oily and combination skin, the procedure can be performed 1–2 times a week.
- For dry facial skin — once a week.

Switching On and Off


- Briefly press the power button to switch on the device.
- Press and hold the power button for 1 second to switch off the device.
- If the device is not used after being switched on, it will automatically switch off after 5 minutes.

Charging the Device

1. Connect the charging cable to the device and a power source.
2. During charging, the indicator flashes. When the device is fully charged, the indicator remains steadily lit.
3. The device cannot be used while charging.

Charging Recommendations

- Use only the original charging cable supplied with the device.
- Before charging, make sure that the Type-C port is dry.
- Charge the device using a DC 5 V / 1 A power source or a USB port with equivalent specifications.
- Using chargers with different specifications may result in slow charging, overheating, or improper operation of the device.

Cleaning the Device

- After completing the procedure, switch off the device.
- Remove any cosmetic product residue from the metal spatula using a soft, damp cloth or a cotton pad moistened with chlorhexidine solution.
- After cleaning, wipe the device dry and allow it to dry completely.
- Store the scrubber in a clean, dry place.

Caution: Do not immerse the device in water or allow water to enter the housing or the charging port. Do not use alcohol, abrasive cleaners, or harsh chemicals for cleaning.

Warranty Terms

Dear Customers!
Please keep these instructions in a safe place in order to use the warranty card. The warranty is valid for one year from the date of purchase (or for 2 or 3 years if an extended warranty has been purchased) in the event of a malfunction under normal operating conditions.

The warranty does not cover the following cases:

1. If the product is damaged due to improper use or lack of proper maintenance.
2. If the product is damaged as a result of disassembly or modification.
3. If parts not supplied by our company have been fitted to or used with the device.
4. Damage caused by force majeure circumstances.

In accordance with the applicable terms, servicing provided outside the warranty is subject to an additional charge.
Any parts of the product replaced during the warranty period remain the property of MEDICA+.

Technical Specifications

Housing material: ABS plastic, stainless steel.
Battery capacity: 450 mAh.
Charging port: USB Type-C.
Power supply: DC 5 V / 1 A.
Overcharge protection: Yes.
Full charging time: approximately 1.5 hours.
Operating time: 45 minutes.
Power: 1.5 W.
Noise level: 65–70 dB.
Weight: 80 g.
Dimensions: 175 × 38 × 14 mm.
Ultrasonic vibration frequency for each mode:
ION−: 1 Hz
ION+: 1 Hz
Lifting: 500 Hz
Sonic: 26 kHz

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