ENG / UA/ PL

PROWATER STATION 8.0

medica+
JAPAN TECHNOLOGY

medica+
JAPAN TECHNOLOGY

ENG

medica+
JAPAN TECHNOLOGY

Instruction

PROWATER STATION 8.0

Contents

Safety Instructions 01-02

Product Description..... 02-04

How to Use..... 05-07

Operating Modes 08-09

Post-Use Care 09-10

Charging 10-11

Daily Maintenance 11

Warranty Terms 12

Technical Specifications 13

Dear Customer,

Thank you for choosing MEDICA+ products. Please read this manual carefully before using the device.

Safety Instructions


1. This model of irrigator should not be used by individuals with reduced physical, sensory, or mental capabilities, unless they are supervised or have been instructed on how to use the device safely.
2. Individuals who have undergone oral or gum surgery within the past two months are advised to consult a dentist before using the device.
3. Do not use the charger outdoors, in damp environments, or near heat sources or hot surfaces. Avoid contact between the charger and water.
4. The manufacturer is not responsible for any damage to the device or harm to the user resulting from misuse or failure to follow the instructions in this manual.
5. Do not use the device for any purpose other than oral cleaning. Misuse may lead to injury, especially if water comes into contact with the eyes, nose, ears, or throat.


6. If you experience pain or discomfort, stop using the device immediately and consult a doctor.
7. Do not fill the water tank with water above 40°C (104°F) to avoid burns.
8. If you have a pacemaker or any other implanted electronic device, consult your doctor before using the product.


Product Description

The MEDICA+ PROWATER STATION 8.0 countertop oral irrigator provides effective and gentle oral cleaning thanks to a powerful pulsating water stream. The device features 4 operating modes with adjustable intensity, allowing you to customize the cleaning process to your individual needs.

A 550 ml water reservoir and a pulsation frequency of 1500–2000 pulses per minute ensure thorough removal of food debris and plaque between teeth and along the gum line, significantly improving overall oral hygiene.

1
2
3

4
5
6


7
8
9
10
11
12
13

1. Power On/Off Button
2. Mode Selection Button
3. Digital Display
4. Irrigator Nozzle
5. Nozzle Release Button
6. Charging Port
7. Water Reservoir
8. Reservoir Lid
9. Nozzle Storage Case
10. Power Switch ON/OFF
11. Handle (Holder)
12. Water Hose
13. Base Unit


We recommend starting with the lowest intensity level by selecting the Child or Pulse mode, as the first uses may cause unusual sensations, including increased gum sensitivity due to the pulsating water stream.
This is a normal reaction and typically subsides once you become accustomed to the device.

How to Use


- To access the nozzles, gently press the storage case — it will open automatically.
- Pull the case upward to remove it.
- After selecting the appropriate nozzle, place the case back and press down until it clicks into place.
- Insert the nozzle into the hole at the top of the handle until you hear a click.
- To replace the nozzle, first turn off the device, then lightly press the release button — the nozzle will detach easily.


The water tank can be filled in two ways:
First method:Disconnect the charging cable, open the reservoir lid while holding the irrigator in a vertical position, pour water into the tank, and securely close the lid (B-2).

B-2

Second method: Remove the tank from the base, open its lid, fill the tank with water, close the lid, and place the tank back onto the base (B-3). Before removing the tank, always ensure the device is powered off. After reinstalling, press the tank down firmly so that no gap remains between it and the base.

B-3

For first-time use, it is recommended to start with the lowest power level and gradually increase intensity over a few weeks. This helps the gums adapt and improves cleaning efficiency. Always place the nozzle in your mouth before turning on the device. Lean over the sink and slightly close your lips to prevent splashing, while allowing water to flow out freely. Switch on the device, and after completing the cleaning process, turn it off (B-4).

B-4

For best results, begin cleaning from the back teeth and move gradually toward the front. Do not direct the water stream straight into periodontal pockets. Aim the stream at a 90° angle to the gum line, gently moving the nozzle along the gums and pausing briefly between each tooth. Continue cleaning until all inner and outer surfaces of the upper and lower teeth are treated.

Operating Modes

The device features four operating modes:

- Strong,
- Normal,
- Pulse
- Child.
 The Strong, Normal, and Child modes support nine intensity levels each. The device automatically remembers the last used intensity level for each of these modes. The Pulse mode includes three intensity levels.

Strong Mode

Designed for intensive cleaning between teeth and removing stubborn plaque. It delivers high water pressure and is suitable for experienced users with healthy gums. Recommended for use after full adaptation to the device.

Normal Mode

An everyday mode that combines efficiency and comfort. It offers an optimal balance between water pressure and gentle action, making it ideal for most users and for regular oral hygiene.

Pulse Mode
Operates as rhythmic water pulses that massage the gums and improve circulation. Recommended for gum care, heightened sensitivity, or during the adjustment period. Also suitable for refreshing cleanings.

Child Mode
A low-intensity mode specially developed for children and users with sensitive gums. Provides gentle and safe cleaning, suitable for beginners and daily delicate oral care.

Post-Use Care

- After use, turn off the device using the power button.
- The irrigator is equipped with an automatic shut-off timer for user safety and device protection.
- In Strong / Normal / Child modes, the device will automatically power off after 4 minutes of continuous operation. In Pulse mode, it powers off after 10 minutes. If the timer activates but you wish to continue the procedure, simply press the power button again.
- The water tank is designed for one cleaning cycle; refill if necessary.
- The memory function stores the last-used mode and intensity for your next session.

- After each use, turn off the device, remove and empty the water tank. Detach the nozzle by pressing the release button on the back of the handle, then rinse it thoroughly under running water. After cleaning, reinsert the nozzle into the handle and wipe the handle and water hose with a dry cloth to remove any remaining moisture.
- To prevent bacterial growth, ensure no water remains inside the nozzle. If the device will not be used for an extended period, store it completely dry.

Charging

- This device is powered by a lithium-ion battery. Before charging, ensure that both the device and your hands are completely dry. Connect the USB charging cable to a power adapter, plug the adapter into an electrical outlet, and then insert the cable's connector into the device's charging port, making sure the connection is secure.
- If the red indicator light turns on after use or the device fails to power on, this indicates a low battery level and that charging is required. A full charging cycle takes approximately 4 hours. While charging, the red indicator light will remain on. Once charging is complete, the green indicator will light up, indicating a full charge.

The device features an automatic power-off function upon reaching full charge to prevent overcharging and protect battery health.

- Use a power adapter with an output of 5V ⎓ 1A during charging. Do not use charging devices with an output exceeding 15W, as this may damage the device or battery.

Daily Maintenance

1. Clean the device before first use and at least once every 6 months if it has not been used for an extended period.
2. Use clean water or a neutral detergent for cleaning. Do not use abrasive, corrosive, or aggressive cleaning agents to avoid damaging the device's surface.
3. Avoid using water with a temperature above 40°C (104°F), as this may cause deformation or damage to the components.
4. After cleaning, thoroughly dry the device and its components before storing.
5. Store the device in a dry, well-ventilated area, away from direct sunlight and heat sources.

Warranty Terms

Dear customers, Please keep this manual in a safe place in order to use it as a warranty certificate.
The product comes with a one-year warranty from the date of purchase (or 2 or 3 years with extended warranty) in case of malfunction under normal usage conditions.

The warranty does not cover the following cases:

- Damage caused by improper use or lack of proper maintenance
- Damage resulting from disassembly or modification
- Use of components or accessories not provided by our company
- Damage caused by force majeure circumstances

According to the warranty terms, any out-of-warranty maintenance is subject to an additional service fee.

Any product parts replaced under the warranty period remain the property of MEDICA+.

Technical Specifications

- Product dimensions: 191.8 × 196.3 × 96.3 mm
- Weight: 926 g
- Pressure range: 80–160 PSI
- Battery capacity: Li-ion, 1500 mAh
- Number of modes: 4 modes
- Auto shut-off: Strong / Normal / Child modes 4 minutes, Pulse mode — 10 minutes
- Intensity levels: Strong / Normal / Child — 9 levels, Pulse — 3 levels
- Input voltage: DC 5 V
- Rated operating power: ≤ 25 W
- Rated input power: ≤ 5 W
- Water tank capacity: 550 ml
- Water pulse frequency: 1500–2000 pulses/min
- Noise level: ≤ 70 dB
- Charging time: 4 hours
- Waterproof rating: IPX7
- Charging adapter: Type-C
- Battery life per charge: up to 15 uses
- Water flow rate: ≥ 290 ml/min

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