ENG / UA/ PL

SkinLifting PRO 4v1

medica+
JAPAN TECHNOLOGY

ENG

medica+
JAPAN TECHNOLOGY

Instruction manual

MEDICA+ SkinLifting PRO 4v1
Lymphatic Drainage Device

Content

Safety Precautions 01-02

Device Description 03

Mode Selection 04-06

Operating Instructions 07-11

Usage Protocol 12-16

Contraindications 17-18

Cleaning the Device 19

Charging and Operating Conditions 20-21

Warranty Terms 22

Technical Specifications 23

Dear Customer!

Thank you for choosing MEDICA+ products. Please read these instructions carefully before using the device.

Safety Precautions

Before use, carefully read these instructions and follow the recommendations provided. Improper use of the device may result in injury, skin irritation, or damage to the product.

- When using the LED mode, light may reach the eye area. Do not look directly at the LEDs. When treating the area around the eyes, it is recommended to keep your eyes closed.
- Do not immerse the device in water or use it while taking a shower or bath.
- Do not charge or use the device near an open flame, high-temperature sources, or in direct sunlight. The built-in lithium battery may overheat or malfunction.
- Keep the device out of the reach of children.
- Use the massager only in accordance with these instructions.
- If you experience pain, discomfort, skin irritation, or any unusual symptoms


during use, stop using the device immediately and consult a doctor.

Charging Safety Precautions

- Use only an AC 100–240 V power source and a DC 5 V / 1 A charger.
- Do not use chargers or adapters that do not meet the specified ratings.
- Do not use a damaged power cable or adapter, or a faulty electrical outlet.
- Do not charge the device in a bathroom or other areas with high humidity.
- Do not touch the charging cable, adapter, or plug with wet hands.
- If a malfunction or unusual situation occurs while the device is charging or operating, disconnect it from the power source immediately.

Device Description

The MEDICA+ SkinLifting PRO 4-in-1 is a multifunctional beauty device designed for comprehensive at-home skincare. It combines advanced technologies such as EMS microcurrents, LED light therapy with blue, red, and green light, vibration massage, and heating and cooling functions. The device is designed to improve skin tone, reduce the appearance of wrinkles, relieve puffiness, and enhance facial skin elasticity.

Mode Selection

Blue Light (460 nm)

Blue light penetrates the upper layers of the skin to a depth of approximately 0.5 mm. It helps soothe the skin and supports the care of problem-prone skin that is susceptible to breakouts. It helps maintain normal sebaceous gland function and is suitable for sensitive skin. It also helps improve the appearance of post-acne skin.

Red Light (630 nm)

Red light penetrates the skin to a depth of 6–8 mm, reaching the dermal layer. It helps stimulate fibroblasts, improve microcirculation and lymphatic flow, support the skin's natural elasticity, even out skin tone, make the skin smoother and more hydrated, and reduce the visible appearance of pigmentation spots.

Green Light (520 nm)

Green light penetrates the upper layers of the skin to a depth of approximately 1–2 mm. It helps even out skin tone, reduce the visible appearance of pigmentation and redness, soothe the skin, and give it a fresher and more radiant appearance. It is suitable for all skin

types, especially sensitive and dull skin.

Electrical Muscle Stimulation (EMS)

EMS technology uses low-level microcurrents to stimulate the facial muscles. The micro-impulses activate neuromuscular tissue, causing gentle muscle contractions. This helps improve skin tone, promote a more defined facial contour, support skin firmness, and reduce the visible appearance of fine lines.

During EMS operation, a mild tingling sensation or muscle contraction may be felt; this is a normal response. If the sensation causes pain or significant discomfort, stop using the device and reduce the intensity level.

Intermittent Vibration (VR)

The vibration massage mode can be used in combination with any other mode. The vibration provides a gentle massage effect, helps improve microcirculation, promotes more even distribution of cosmetic products, and increases comfort during EMS use.

Cooling Mode


The cooling temperature is 15 ± 3 °C. The

cooling effect helps tone and firm the skin, tighten pores and small capillaries, and reduce puffiness. Cooling also helps retain moisture in the skin and reduces the likelihood of impurities re-entering the pores.

Heating Mode

The heating temperature is 42 ± 3 °C. Gentle heat stimulates microcirculation and metabolic processes in the skin, promotes more even distribution of cosmetic products across the skin's surface, and improves their absorption by the stratum corneum.

Operating Instructions

Button 2: Switching the Device On/Off and Selecting LED Light Therapy Modes

Switching on the device

Press and hold Button 2 for 1–2 seconds. After the audible signal, the device will switch on and enter standby mode, and all indicators will light up.

Note: If no action is performed for 1 minute, the device will emit two audible signals and switch off automatically.

Switching off the device

Press and hold Button 2 for 1–2 seconds. After the audible signal, the device will switch off.

Selecting LED light therapy modes

After switching on the device, briefly press Button 2 to select the desired light therapy mode:

- 1 press – blue light (Blue Light);
- 2 presses – red light (Red Light);
- 3 presses – green light (Green Light);
- 4 presses – LED light therapy off.

Each light therapy mode operates for 10 minutes, after which it ends automatically.

Warning: Do not use LED light therapy if you have photosensitivity or are taking medications that increase sensitivity to light.

Button 1: EMS Intensity Selection and Vibration Mode Control

Adjusting the EMS Intensity

After switching on the device, briefly press Button 1 to activate the EMS (electrical muscle stimulation) function. Intensity Level 1 is activated by default, and the corresponding indicator lights up.

Each subsequent brief press of Button 1 switches the intensity level:

- 1 press – Level 1;
- 2 presses – Level 2;
- 3 presses – Level 3.

Note: Once activated, the EMS function cannot be switched off separately. To deactivate it, switch the device off and then switch it on again.

Gear 1

Gear 2

Gear 3

Switching the Vibration Mode (VR) On and Off

- Press and hold Button 1 to switch on the vibration mode.
- Press and hold Button 1 again to switch off the vibration mode.

The vibration mode can be used separately or in combination with any other operating mode of the device.

Note: The continuous operating time of the EMS and vibration (VR) functions is 10 minutes, after which they switch off automatically.

Button 3: Switching Between Heating and Cooling Modes

The cooling or heating mode can operate simultaneously with the EMS, vibration (VR), and LED light therapy functions.

Note: The cooling and heating modes cannot operate simultaneously.

Cooling Mode

1 brief press of Button 3 – switches on the cooling mode. The operating time is 3 minutes.

Heating Mode

2 brief presses of Button 3 – switch on the heating mode. The operating time is 5 minutes.

3 brief presses of Button 3 – switch off the cooling and heating modes.

After the activated mode has finished operating, the device enters standby mode.

When switching between modes, the timer for the selected mode restarts.

If the heating or cooling mode has ended while other functions—EMS, vibration, or LED light therapy—are still operating, the device will enter standby mode only after all active functions have finished.

Usage Protocol

Before using the device for the first time, always perform a preliminary test to ensure that your skin tolerates the procedure well. This will help prevent possible irritation or adverse reactions.

Recommended Treatment Frequency

Intensive course: 3–4 treatments per week.
Maintenance course: After achieving the desired result, it is recommended to use the device 1–2 times per week to maintain the effect.

1. Before each use, thoroughly cleanse the facial skin of makeup and impurities. After cleansing, apply a suitable gel, booster, or serum to the skin.
2. When using the EMS (electrical muscle stimulation) mode, it is recommended to apply MEDICA+ Conductive Lifting Gel or another water-based product beforehand. This improves microcurrent conductivity and ensures more effective operation of the device.

Mode Combinations

The protocols below are basic

recommendations. Depending on your skin's needs, you can combine the device's operating modes and LED light therapy modes to create an individual skincare protocol.

Basic Care

Maintaining skin cleanliness, reducing puffiness, and providing a mild toning effect.
Mode combination:

- Blue LED light
- VR (vibration) — optional
- Cooling mode

How to Use:

1. Cleanse the skin and apply a water-based booster or serum.
2. Switch on the blue LED light.
3. Activate VR for a gentle massage.
4. Add the cooling mode for 3 minutes.

This combination is used to soothe the skin, reduce inflammation, tighten pores, and decrease puffiness.

Anti-Aging Care

Improving skin firmness, defining the facial contour, and reducing the appearance of wrinkles.
Mode Combination:

- EMS, Level 1–3
- Red LED light
- VR vibration
- Heating mode

How to Use:

1. Apply a conductive gel.
2. Activate EMS and select a comfortable intensity level.
3. Switch on the red LED light.
4. Add the heating mode for 5 minutes.
5. Activate VR as needed for an additional massage effect.

This combination helps improve skin tone and firmness, reduce the visible appearance of wrinkles, and stimulate microcirculation.

Deep Hydration

Deep hydration and enhanced effectiveness of skincare products.

Mode Combination:

- EMS at a low intensity level
- Red LED light
- Heating mode

How to Use:

1. Apply a serum or a special booster for light therapy.
2. Activate EMS at the minimum intensity level.

3. Switch on the red LED light.
4. Add the heating mode for 5 minutes.

This combination is used to improve the absorption of cosmetic products, provide deep hydration, and enhance their effectiveness.

Evening Out Skin Tone

Helps even out skin tone, reduce the visible appearance of pigmentation spots, and improve the skin's radiance and freshness.

Mode Combination:

- Green LED light
- EMS at a low intensity level
- Heating mode
- VR vibration — optional

How to Use:

1. Apply a serum or a special booster for light therapy.
2. Switch on the green LED light.
3. Activate EMS at the minimum intensity level.
4. Add the heating mode for 5 minutes.
5. Activate VR as needed for an additional massage effect.

This combination helps even out skin tone, enhance the skin's natural radiance, reduce the visible appearance of pigmentation, and maintain a more even complexion.

During the procedure, move the device smoothly along the massage lines without holding it in one area for an extended period. The duration of one treatment should not exceed 10 minutes.

Contraindications

Do not use the device if you:

- Have a pacemaker or other implanted electronic medical devices;
- Use life-support equipment;
- Have diseases of the internal organs in an acute stage within the treatment area;
- Have facial nerve disorders;
- Have mental health disorders;
- Are in the acute stage of any illness;
- Have uncontrolled or abnormally high or low blood pressure;
- Have an infectious disease;
- Have tuberculosis;
- Have a blood disorder;
- Have a malignant tumour;
- Are pregnant, breastfeeding, have a fever, or are undergoing dental treatment, unless approved by a doctor;
- Have atopic dermatitis, other atopic skin conditions, or a tendency to allergic reactions;
- Have highly sensitive skin;
- Have been taking steroid hormone medications for a prolonged period or have vascular skin changes associated with liver disease;

- Have epilepsy;
- Have a heart rhythm disorder;
- Have recently received botulinum toxin injections or dermal fillers, until approved by a doctor;
- Have recently undergone aggressive cosmetic procedures, until the skin has fully recovered.

Do not use the device on the following areas:

- Open wounds;
- Areas that have undergone plastic surgery until fully healed;
- Directly on the eyeballs or upper eyelids;
- Inside the mouth;
- Mucous membranes;
- The intimate area;
- Areas with severe burns;
- Areas with redness, itching, swelling, or inflammation;
- Areas with pain or impaired sensitivity;
- Areas where metal, silicone, or plastic implants are located beneath the skin.

Cleaning the Device

1. Before cleaning the device, always switch it off and disconnect it from the power source.
2. Clean the massage head and the device body using a damp antiseptic wipe or a cotton pad moistened with chlorhexidine. Note: The device has an IPX4 protection rating. The massage head is water-resistant; however, the device itself must not be immersed in water or rinsed under running water.
3. After cleaning, store the device in a dry, cool place protected from direct sunlight.
4. Do not use alkaline products, alcohol, volatile oils, solvents, or other aggressive cleaning agents.
5. Clean the device regularly and promptly after each use.
6. If the device will not be used for an extended period, be sure to clean it before storage.
7. Do not charge the device while it is wet.
8. Do not use oily or abrasive cosmetic products during treatments.

Charging and Operating Conditions

1. The device is not intended for use while charging. Do not switch it on while it is connected to a power source.
2. Before first use, it is recommended to fully charge the device to ensure optimal performance and a long battery service life.
3. To maintain stable battery performance, charge the device at least once a month. Failure to follow this recommendation may result in deep discharge and damage to the battery.
4. Low battery indication. When the battery level is low, the device provides the following alerts:
 - The indicator flashes rapidly 10 times;
 - The audible signal sounds 10 times;
 - The device then switches off automatically.
5. It is recommended to charge the device at a room temperature of 10°C–30°C.
6. The device can be charged in one of two ways:
 - Connect the charging cable directly to the device and the power adapter.

- Place the device on the charging base and connect the cable to the base and the power adapter.
 Note: The device operates only from a DC 5 V power source.

7. During charging, the indicator is activated. Once the device is fully charged, the indicators remain continuously illuminated.
8. Do not charge the device in a bathroom or under high-humidity conditions, as this may result in electric shock or a short circuit.

Warranty Terms

Dear Customers,
Please keep this instruction manual in a safe place so that you can use the warranty card. The warranty is valid for one year from the date of purchase, or for two or three years if an extended warranty has been purchased, in the event of a malfunction under normal operating conditions.
The warranty does not cover the following cases:

1. If the product is damaged due to improper use or lack of proper maintenance.
2. If the product is damaged as a result of disassembly or modification.
3. If parts not supplied by our company have been installed or used with the device.
4. Damage caused by force majeure circumstances.

In accordance with the applicable terms, servicing provided outside the warranty period is subject to an additional charge.
Parts replaced during warranty repair become the property of the MEDICA+ service center.

Technical Specifications

Housing material: PC plastic, ABS plastic, zinc alloy, silicone
Device dimensions: 179.5 × 72 × 48.4 mm
Charging base dimensions: 81 × 69.8 × 64.8 mm
Device weight: 164 g
Charging base weight: 115 g
Rated power: 10 W. Operating voltage: 3.7 V
Power supply: built-in lithium-polymer battery
Battery capacity: 1350 mAh
Charging time: approximately 3 hours
Battery life: up to 35 minutes
Protection rating: IPX4 (protection against water splashes)
EMS (electrical muscle stimulation): 3 intensity levels

- Level 1 — 0.04 mA
- Level 2 — 0.05 mA
- Level 3 — 0.11 mA

LED light therapy: red light — 630 nm, blue light — 460 nm, green light — 520 nm
Light intensity: red light — 10 mW/cm², blue light — 24 mW/cm², green light — 21 mW/cm²
Vibration speed: 7,000 rpm
Heating: 42 ± 3 °C. Cooling: 15 ± 3 °C
Automatic shut-off: after 10 minutes
Heating mode: 5 minutes
Cooling mode: 3 minutes

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